


Umps Pendant Chain

Attach your pendant

Remove the existing lanyard first. If you need help removing it or taking your pendant out of a wrist case, ask your village coordinator, carer or Umps.

1 Open the clip

Press and hold the small lever on the pendant attachment clip to open it.



2 Attach and close

Pass the open clip through the eyelet at the top of your pendant. Release the lever and check the clip closes fully around the eyelet.



The magnetic clasp is a safety breakaway. Do not glue, tie or otherwise prevent it from opening.

Cleaning and care

After showering, pat the chain and clasp dry with a soft towel. Clean with a soft cloth, lukewarm water and mild soap if needed. Rinse briefly and pat dry. Avoid soaking, abrasive cleaners, silver polish and jewellery-cleaning machines. The silver finish may darken over time. Do not scrub it to restore shine.

Check for wear

Check the chain, clips, magnetic clasp and pendant eyelet occasionally. Contact Umps if anything is worn, rough or damaged, or does not close securely.

The chain contains nickel. If your skin becomes itchy, red or sore, stop using the chain and contact Umps for another wearing option.

Need support? Contact Umps.

 1300 907 272

 support@umps.com

 www.umps.com

View our complete user manuals, datasheets and product specifications at www.umps.com/how-it-works